








PlastPharm packaging catalogue







VIALS						
Picture	Description	Neck finish	Height	OD	Material Glass	Color
	VIAL 2 ml	ISO 13	35 mm	16 mm	Type I	Clear
	VIAL 6 ml	ISO 20	40 mm	22 mm	Type I	Clear / Amber
	VIAL 7ml	ISO 20	40.75 mm	22 mm	Type I ; II ; III	Clear / Amber





PlastPharm packaging Catalogue

	<p style="text-align: center;">VIAL 7 ml</p>	<p style="text-align: center;">ISO 20</p>	<p style="text-align: center;">38 mm</p>	<p style="text-align: center;">21 mm</p>	<p style="text-align: center;">Type III</p>	<p style="text-align: center;">Clear</p>
	<p style="text-align: center;">VIAL 8 ml</p>	<p style="text-align: center;">ISO 20</p>	<p style="text-align: center;">46.8 mm</p>	<p style="text-align: center;">23 mm</p>	<p style="text-align: center;">Type I ; II ; III</p>	<p style="text-align: center;">Clear</p>
	<p style="text-align: center;">VIAL 8 ml</p>	<p style="text-align: center;">ISO 20</p>	<p style="text-align: center;">45 mm</p>	<p style="text-align: center;">22 mm</p>	<p style="text-align: center;">Type I</p>	<p style="text-align: center;">Clear / Amber</p>
	<p style="text-align: center;">VIAL 10ml QUARTER</p>	<p style="text-align: center;">ISO 20</p>	<p style="text-align: center;">49 mm</p>	<p style="text-align: center;">24.5 mm</p>	<p style="text-align: center;">Type III</p>	<p style="text-align: center;">Amber</p>





PlastPharm packaging Catalogue

	<p>VIAL 10ml</p>	<p>ISO 20</p>	<p>53.5 mm</p>	<p>25.4 mm</p>	<p>Type I ; II ; III</p>	<p>Clear / Amber</p>
	<p>VIAL 10ml</p>	<p>ISO 20</p>	<p>49 mm</p>	<p>24.5 mm</p>	<p>Type I ; III</p>	<p>Clear / Amber</p>
	<p>VIAL 10ml</p>	<p>ISO 20</p>	<p>45 mm</p>	<p>24 mm</p>	<p>Type I</p>	<p>Clear / Amber</p>
	<p>VIAL 10ml</p>	<p>ISO 20</p>	<p>45 mm</p>	<p>24 mm</p>	<p>Type I</p>	<p>Clear / Amber</p>





PlastPharm packaging Catalogue

	<p style="text-align: center;">VIAL 20 ml</p>	<p style="text-align: center;">ISO 20</p>	<p style="text-align: center;">28.8 mm</p>	<p style="text-align: center;">26.5 mm</p>	<p style="text-align: center;">Type I ; III</p>	<p style="text-align: center;">Clear / Amber</p>
	<p style="text-align: center;">VIAL 15 ml</p>	<p style="text-align: center;">ISO 20</p>	<p style="text-align: center;">62.5 mm</p>	<p style="text-align: center;">24.8 mm</p>	<p style="text-align: center;">Type III</p>	<p style="text-align: center;">Amber</p>
	<p style="text-align: center;">VIAL 15 ml</p>	<p style="text-align: center;">ISO 20</p>	<p style="text-align: center;">25 mm</p>	<p style="text-align: center;">32 mm</p>	<p style="text-align: center;">Type I ; II ; III</p>	<p style="text-align: center;">Clear / Amber</p>
	<p style="text-align: center;">VIAL 20 ml</p>	<p style="text-align: center;">ISO 20</p>	<p style="text-align: center;">55 mm</p>	<p style="text-align: center;">30 mm</p>	<p style="text-align: center;">Type I</p>	<p style="text-align: center;">Clear / Amber</p>





PlastPharm packaging Catalogue

	VIAL 15 ml	ISO 20	52 mm	26.5 mm	Type III	Clear
	VIAL 20 ml	ISO 20	55 mm	30 mm	Type I	Clear / Amber
	VIAL 25 ml	ISO 20	65 mm	30 mm	Type III	Clear
	VIAL 30 ml	ISO 20	67.1 mm	31.7 mm	Type I	Clear / Amber




PlastPharm packaging Catalogue

	VIAL 30 ml	ISO 20	65.5 mm	32 mm	Type I	Clear / Amber
	VIAL 30 ml	ISO 20	75 mm	30 mm	Type I	Clear
	VIAL 30 ml	ISO 20	62.8 mm	36 mm	Type I	Clear / Amber
	VIAL 45ml	ISO 20	73 mm	36.3 mm	Type I	Clear / Amber

PlastPharm packaging Catalogue

	<p>Infusion IV 50ml max level mark 50ml</p>	<p>ISO 32</p>	<p>68 mm</p>	<p>46 mm</p>	<p>Type I ; II</p>	<p>Clear</p>
	<p>VIAL 50 ml</p>	<p>ISO 20</p>	<p>73 mm</p>	<p>42.5 mm</p>	<p>Type I</p>	<p>Clear / Amber</p>
	<p>Infusion IV 60 ml</p>	<p>ISO 32</p>	<p>76.2 mm</p>	<p>46 mm</p>	<p>Type III</p>	<p>Clear</p>
	<p>VIAL 100 ml</p>	<p>ISO 20</p>	<p>94.5 mm</p>	<p>51.69 mm</p>	<p>Type I</p>	<p>Clear / Amber</p>

PlastPharm packaging Catalogue

	<p>Infusion IV 100ml max level mark 100ml</p>	<p>ISO 32</p>	<p>104 mm</p>	<p>49 mm</p>	<p>Type I ; II ; III</p>	<p>Clear</p>
	<p>Infusion IV 125ml max level mark 125ml</p>	<p>ISO 32</p>	<p>98 mm</p>	<p>54.4 mm</p>	<p>Type III</p>	<p>Clear</p>
	<p>Infusion IV 200ml mark 200ml</p>	<p>ISO 32</p>	<p>118 mm</p>	<p>66 mm</p>	<p>Type I</p>	<p>Clear</p>